

90D Aspheric Lens

This slit lamp lens is used for general examination of the fundus and small pupil examinations.



Product Description

Product Code	Description	Field of View
MIPL/K3	90D Aspheric Lens	94°

Features

- This Slit lamp lens is the most popular power for non-contact fundus examination.
- It offers a wide-field fundus imaging out to and beyond the equator.
- A high magnification provides an increased illuminated field area which makes it ideal for posterior pole scanning.
- Made in glass.
- Aspheric Optics
- CE Certified.

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase-1, New Delhi -110020, India

T: (+91-11) 4270 1030, 31 | E: Sales@madhuinstruments.com | W: www.madhuinstruments.com

Customer Care No. +91 9810866176