




Irrigating Vitrectomy Lens Set

This set comprises of 3 lenses which have a provision to attach the infusion line for irrigation during a VR surgery.



Product Description

Product Code	Description	Field of View	Qty. per box
MIPL/L2/BC	-90D/Biconcave		1 pc.
MIPL/L2/WF	Widefield		1 pc.
MIPL/L2/CF	Central Fundus		1 pc.

Features

- High quality glass optics for clear images
- Lens description engraved on the handle for easy identification
- Flexibility to order individual lenses
- CE Certified

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase-1, New Delhi -110020, India

T: (+91-11) 4270 1030, 31 | E: Sales@madhuinstruments.com | W: www.madhuinstruments.com

Customer Care No. +91 9810866176