



Lens Glide

Disposable lens glide is used for insertion of IOLs or other microsurgical instruments in the anterior chamber.



Product Description

Product Code	Description	Qty. per box
MIPL/D5	Lens Glide	10 Pcs. (Sterile)

Features

- Lens glide is a double ended Polypropylene device
- The narrow end is used during phacoemulsification surgical procedure whereas the broader end is used during non-phaco surgical procedures.
- CE certified

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase-1, New Delhi -110020, India

T: (+91-11) 4270 1030, 31 | E: Sales@madhuinstruments.com | W: www.madhuinstruments.com

Customer Care No. +91 9810866176

Ref. No.: MIPL/Pub/SPB/D5/01, Issue No.: 01, Issue Date: 01.04.2020, Rev. No. 01, Rev. Date: 01.09.2020