



## Peroperative Keratoscope

Peroperative keratoscope is used during an ophthalmic surgical procedure to judge the astigmatism induced in the cornea while putting sutures.



### Product Description

Product Code	Description	Qty. per pack
MIPL/N5	Peroperative Keratoscope	1 Pc. (Unsterile)

### Features

- This is a mechanical device made out of anodized aluminum tube.
- The inner wall of the tube has alternate black and white rings of different thickness.
- The reflection of these rings on the patient's cornea helps in analyzing the astigmatism induced due to tight and loose sutures.
- Can be resterilised by autoclave or ETO
- CE certified

## Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

## Cleaning and Disinfection

For cleaning and disinfection of the device, kindly refer to the IFU supplied with the device or contact our customer care number



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase-1, New Delhi -110020, India

T: (+91-11) 4270 1030, 31 | E: Sales@madhuinstruments.com | W: www.madhuinstruments.com

Customer Care No. +91 9810866176