



PVA Instrument Wipe

High quality PVA sponge wipe for care and cleaning of extremely delicate micro instruments.



Product Description

Product Code	Description	Qty. per box
MIPL/B5	PVA Instrument Wipe	1 pc. (Sterile)

Features

- Special PVA material is highly absorbent and particulate free
- Does not damage extremely delicate instruments
- Prolongs the life of these expensive instruments

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through gamma radiation according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)