

78D Aspheric Lens

This slit lamp lens is used for general examination of the fundus.



Product Description

Product Code	Description	Field of View	Qty. per box
MIPL/K2	78D Aspheric Lens	84°	1 pc.

Features

- This Slit lamp lens offers an ideal combination magnification and field of view.
- It also provides an extended working distance from cornea.
- Superior glass optics offers high clarity images.
- Made in glass.
- Aspheric Optics
- CE Certified.

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

Cleaning & Disinfection

For cleaning and disinfection of the lens, kindly refer to the IFU supplied with the lens or contact our customer care.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)