


## Single Mirror Goniolens

This gonio lens is used for viewing the anterior chamber angle and central retina.



### Product Description

Product Code	Description	Field of View	Qty. per box
MIPL/I1	Single Mirror Goniolens	62°	1 pc.

### Features

- Due to its small size, a single mirror goniolens is very convenient to handle.
- The mirror angle is 62° which when rotated through 360° gives a clear view of the anterior chamber angle.
- Its small contact surface makes it very useful while diagnosing adults as well as children with small palpebral fissure.
- Made in glass.
- Contact medium required.
- CE Certified.

### Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

### Cleaning & Disinfection

For cleaning and disinfection of the lens, kindly refer to the IFU supplied with the lens or contact our customer care



**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)