

90D Aspheric Lens

This slit lamp lens is used for general examination of the fundus and small pupil examinations.



Product Description

Product Code	Description	Field of View	Qty. per box
MIPL/K3	90D Aspheric Lens	94°	1 pc.

Features

- This Slit lamp lens is the most popular power for non-contact fundus examination.
- It offers a wide-field fundus imaging out to and beyond the equator.
- A high magnification provides an increased illuminated field area which makes it ideal for posterior pole scanning.
- Made in glass.
- Aspheric Optics

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)