


## Silicone Tip Cannula

A disposable silicone tip cannula is used for atraumatic removal of intraocular fluids by active or passive aspiration during a VR surgery.



### Product Description

Product Code	Description	Silicon Tip Length	Qty. per box
MIPL/E2/STR/20G	Silicone Tip Cannula (Straight/Regular) – 20G	1 mm	10 pcs (Sterile)
MIPL/E2/STBR/20G	Silicone Tip Cannula (Straight/Brush) – 20G	1 mm	10 pcs (Sterile)
MIPL/E2/STR/23G	Silicone Tip Cannula (Straight/Regular) – 23G	1 mm	10 pcs (Sterile)
MIPL/E2/STBR/23G	Silicone Tip Cannula (Straight/Brush) – 23G	1 mm	10 pcs (Sterile)
MIPL/E2/STR/25G	Silicone Tip Cannula (Straight/Regular) – 25G	1 mm	10 pcs (Sterile)
MIPL/E2/STBR/25G	Silicone Tip Cannula (Straight/Brush) – 25G	1 mm	10 pcs (Sterile)
MIPL/E2/GN/STR/23G	Silicone Tip Cannula (Straight/Regular) Green - 23G	1 mm	10 pcs (Sterile)
MIPL/E2/GN/STR/25G	Silicone Tip Cannula (Straight/Regular) Green - 25G	1 mm	10 pcs (Sterile)

- Ultra thin wall cannula for maximum flow
- CE Certified

Product Code	Description	Silicon Tip Length	Qty. per box
MIPL/E2/BTC/20G	Blunt Tip Cannula – 20G	N/A	10 pcs (Sterile)
MIPL/E2/BTC/23G	Blunt Tip Cannula – 23G	N/A	10 pcs (Sterile)
MIPL/E2/BTC/25G	Blunt Tip Cannula – 25G	N/A	10 pcs (Sterile)

- Ultra thin wall cannula for maximum flow



## Features

- Ultra thin wall cannula for maximum flow
- CE certified

## Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

## Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: [sales@madhuinstruments.com](mailto:sales@madhuinstruments.com) | [www.madhuinstruments.com](http://www.madhuinstruments.com)

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)