

Lacrimal Cannula

Disposable ophthalmic cannula used for treatment of canaliculus system stentosis or nasolacrimal duct obstruction.



Product Description

Product Code	Description	Qty. per box
MIPL/E5/23R/7.5	Lacrimal Intubation Set - 7.50cm X 23G (Regular)	10 pcs. (Sterile)
MIPL/E5/23O/7.5	Lacrimal Intubation Set - 7.50cm X 23G (Olive)	10 pcs. (Sterile)
MIPL/E5/23R/11	Lacrimal Intubation Set - 11.00cm X 23G (Regular)	10 pcs. (Sterile)
MIPL/E5/23O/11	Lacrimal Intubation Set - 11.00cm X 23G (Olive)	10 pcs. (Sterile)
MIPL/E5/23R/15	Lacrimal Intubation Set - 15.00cm X 23G (Regular)	10 pcs. (Sterile)
MIPL/E5/23O/15	Lacrimal Intubation Set - 15.00cm X 23G (Olive)	10 pcs. (Sterile)
MIPL/E5/23R/17.5	Lacrimal Intubation Set - 17.50cm X 23G (Regular)	10 pcs. (Sterile)
MIPL/E5/23O/17.5	Lacrimal Intubation Set - 17.50cm X 23G (Olive)	10 pcs. (Sterile)
MIPL/E5/20STR/4.5	Lacrimal Intubation Set - DCR Straight - 4.50cm x 20G	10 pcs. (Sterile)
MIPL/E5/20ANG/4.5	Lacrimal Intubation Set - DCR Angled - 4.50cm x 20G	10 pcs. (Sterile)

Features

- High quality stainless steel probes and silicon tubes used
- EO Sterile

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)