


Hydrodissection/Hydrodelineation Cannula

Disposable ophthalmic cannula used for effectively separate layers of nucleus and cortex in cataract surgery procedures.



Product Description

Product Code	Descripton with Gauge	Needle Length	Bend Length	Qty. per box
MIPL/E9/A3	Hydrodissection/Hydrodelineation Cannula - 23G	22 mm	9 mm	10 Pcs. (Sterile)
MIPL/E9/A1	Hydrodissection/Hydrodelineation Cannula - 25G	22 mm	9 mm	10 Pcs. (Sterile)
MIPL/E9/A4	Hydrodissection/Hydrodelineation Cannula - 27G	22 mm	9 mm	10 Pcs. (Sterile)
MIPL/E9/A2	Hydrodissection/Hydrodelineation Cannula (J-Shaped) - 25G	22 mm	9 mm	10 Pcs. (Sterile)

Features

- Thin wall cannula used for maximum flow
- CE certified

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)