


Corneal Trephine

Nanoedge corneal trephines are ready-to-use circular blades designed for accurate cutting of corneal donor buttons to be used in penetrating keratoplasty.



Product Description

Product Code	Description	Qty. per box
MIPL/C1/2.00	2.00 mm Diameter	1 pc. (Sterile)
MIPL/C1/3.00	3.00 mm Diameter	1 pc. (Sterile)
MIPL/C1/5.00	5.00 mm Diameter	1 pc. (Sterile)
MIPL/C1/6.00	6.00 mm Diameter	1 pc. (Sterile)
MIPL/C1/6.50	6.50 mm Diameter	1 pc. (Sterile)
MIPL/C1/7.00	7.00 mm Diameter	1 pc. (Sterile)
MIPL/C1/7.25	7.25 mm Diameter	1 pc. (Sterile)
MIPL/C1/7.50	7.50 mm Diameter	1 pc. (Sterile)
MIPL/C1/7.75	7.75 mm Diameter	1 pc. (Sterile)
MIPL/C1/8.00	8.00 mm Diameter	1 pc. (Sterile)
MIPL/C1/8.25	8.25 mm Diameter	1 pc. (Sterile)
MIPL/C1/8.50	8.50 mm Diameter	1 pc. (Sterile)
MIPL/C1/8.75	8.75 mm Diameter	1 pc. (Sterile)
MIPL/C1/9.00	9.00 mm Diameter	1 pc. (Sterile)
MIPL/C1/9.50	9.50 mm Diameter	1 pc. (Sterile)
MIPL/C1/10.00	10.00 mm Diameter	1 pc. (Sterile)
MIPL/C1/10.50	10.50 mm Diameter	1 pc. (Sterile)
MIPL/C1/11.00	11.00 mm Diameter	1 pc. (Sterile)
MIPL/C1/11.50	11.50 mm Diameter	1 pc. (Sterile)
MIPL/C1/12.00	12.00 mm Diameter	1 pc. (Sterile)
MIPL/C1/14.00	14.00 mm Diameter	1 pc. (Sterile)
MIPL/C1/16.00	16.00 mm Diameter	1 pc. (Sterile)
MIPL/C1/18.00	18.00 mm Diameter	1 pc. (Sterile)

Features

- An ultrasharp thin profile cutting edge minimizes damage to the endothelial cells
- Designed for a precise and symmetrical cut of the corneal tissue
- See through design helps in accurate positioning
- Can be used with or without an endothelial punch
- Innovative protective packaging to ensure safety of the cutting edge
- EO Sterile
- CE certified

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)