

Eye Shield

Eye shield is used for a short duration to protect the patient's eye from accidental injury or rubbing after a cataract surgery.



Product Description

Product Code	Description	Qty. per box
MIPL/O3/01	Clear polycarbonate Eye shield	100 Pcs. (Unsterile)
	Clear polycarbonate Eye shield	100 Pcs. (Sterile)
MIPL/O3/02	Polypropylene Eye shield	100 Pcs. (Unsterile)

Features

- Made from high quality clear polycarbonate or polypropylene
- The eye shield covers the patient's eye using an adhesive tape or a patch
- Smooth edges for maximum comfort
- CE certified (Sterile only)

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram from where appropriate. The international symbols used on the packaging are also explained.

Sterilization

Eye Shields are supplied both non sterile and sterile as per buyers requirement.



MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)