


Silicon Band

Scleral buckling devices are used in a retinal detachment surgery.





Product Description

Product Code	Description	Qty. per box
MIPL/H2/2.5	2.5 mm (Style 240) 	10 pcs. (Sterile)
MIPL/H2/3.5	3.5 mm (Style 41) 	10 pcs. (Sterile)
MIPL/H2/4.0	4.0 mm (Style 42) 	10 pcs. (Sterile)

SilicoCare

Features

- Made from medical grade silicon in USA
- Double pouch packed
- EO Sterilised

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate . The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.

madhu
INNOVATING POSSIBILITIES

MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)

Ref. No.: MIPL/SPB/H2/01, Rev. No. 04, Rev. Date: 16.09.2026