Opticlear
Ophthalmic Lenses

Vitrectomy Infusion Cannula

Disposable cannula used for infusion of irrigation fluids or surgical fluids like silicon oil/PFCL during vitreo-retina ophthalmic surgical procedures.





4.0 mm



6.0 mm

Product Description

Product Code	Description	Qty. per box
MIPL/E1/4.0	4.0 mm tip (20G)	10 pcs. (Sterile)
MIPL/E1/6.0	6.0 mm tip (20G)	10 pcs. (Sterile)

Features

- Thin wall stainless steel tube is used for maximum flow
- Designed especially for silicon oil injection
- Superior grip of silicon tube at cannula end as well as connector end

Opticlear
Ophthalmic Lenses

Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.

madhu
INNOVATING POSSIBILITIES

MADHU INSTRUMENTS PVT LTD

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)

Ref. No.: MIPL/SPB/E1/01, Rev. No. 04, Rev. Date: 25.09.2026