# Steady Ring ( Self Stabilizing Ring )

A disposable vitrectomy ring to stabilize the vitrectomy lenses without sutures during a VR surgery.






## Product Description

| Product Code | Description | Qty. per box |
|---|---|---|
| MIPL/D13/SI | Steady Ring (Silicone) | 1 Pc. |

## Features

- Special 'tripod design' for better stability during vitreo retinal surgical procedures
- Eliminates the need for a skilled assistant
- No need to suture the lens to the sclera during the surgery
- Supplied Unsterile

## Instructions for Use

The instructions for use and surgical technique are supplied in the form of a leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are also explained.

## Sterilization

Products are sterilized through Ethylene Oxide (EO) according to a validated cycle.

madhu
INNOVATING POSSIBILITIES

**MADHU INSTRUMENTS PVT LTD**

A-260, Okhla Industrial Area, Phase - 1, New Delhi - 110020, India

T: (+91-11) 4270 1030, 31 | E: sales@madhuinstruments.com | www.madhuinstruments.com

Customer Care No. +91 9810866176 (for international queries), +91 7428186111 (for domestic queries)

Ref. No.: MIPL/SPB/D13/01, Rev. No. 02, Rev. Date: 30.01.2025